



The **power** to achieve

Power to achieve

**A Professional Development Program
for the Mountbatten Braille**

Workshop One

Mountbatten Easy Guide

© 2010 - 2012

Harpo Sp. z o. o.

Poznan, Poland

www.mountbattenbraille.com

Humanware Sales & Support

800-722-3393 (US)

819-471-4818 (Canada)

www.humanware.com

Contents

Introduction.....	<u>1</u>
Mountbatten Orientation - Parts.....	<u>2</u>
Location of Parts.....	<u>3</u>
Mountbatten Orientation - Ports.....	<u>3</u>
Location of Ports.....	<u>5</u>
Features for loading paper.....	<u>6</u>
Loading Paper.....	<u>7</u>
Start-Up.....	<u>8</u>
Battery Care.....	<u>9</u>
The Mountbatten Keyboard.....	<u>10</u>
Using the Mountbatten Keyboard.....	<u>11</u>
Turning Commands On and Off.....	<u>12</u>
Why do we use commands?.....	<u>13</u>
Speech Commands	<u>14</u>
Formatting Commands.....	<u>15</u>
File Commands	<u>16</u>
Ways to adjust noise level.....	<u>17</u>
Adjusting the Braille Impact.....	<u>17</u>
Connecting Mountbatten to a PC with MB-Comm.....	<u>18</u>
Chat Mode using the PC with MB-Comm.....	<u>20</u>
Embossing from the PC with MB-Comm.....	<u>21</u>
Embossing Commands	<u>22</u>
Using the Mimic	<u>23</u>
Forward Translation Commands	<u>24</u>
Changing the Braille Code.....	<u>25</u>
The Options Menu	<u>26</u>
The Options Menu - Calculator.....	<u>27</u>
The Options Menu - Music.....	<u>29</u>
Braille Music Notes.....	<u>31</u>
Braille Graphics.....	<u>31</u>
Writing Numbers in a Command	<u>32</u>
Braille Writing Commands	<u>32</u>
Miscellaneous Commands.....	<u>33</u>
Braille characters and contractions	<u>34</u>
Student Evaluation Tool.....	<u>35</u>

Introduction

This guide has been prepared to provide you with a quick reference tool to use:

- during Workshop One
- as a resource when using your Mountbatten

During the Workshop, we suggest you mark the Easy Guide to help you remember those sections covered.

Some information in the Easy Guide is introduced, but not applied. This allows those already familiar with the Mountbatten to access information when needed, and for those beginners who want to, to go ahead!

Suggestions:

Store the Easy Guide near the Mountbatten.

Introduce the Guide and how to use it to all those who may need to access the information, such as the Aide, Class Teacher and the student.

Copy any procedures used regularly and put them on the wall near the Mountbatten for easy reference.

Braille out reference information as the student user learns about the Mountbatten parts, ports and functionality so that they build skills of independence!

Mountbatten Orientation - Parts

Keyboard

Ergonomically designed, yellow keyboard with six black Braille keys and seven blue keys that tell the Mountbatten what to do.

Embossing Head

Silver metal with yellow marker to help indicate position on the page. On start-up, the embossing head moves to the right hand margin and back to identify paper size.

Reading Table

A flat silver metal section on which to press when reading. The table has a groove along it to align the top of the page. Lay a new page flat on the table to assist with loading paper.

Margin Adjusters

Both margin adjusters have grip wheels to help feed the paper through. The left margin adjuster is fixed. The right margin adjuster moves to accommodate paper size.

Paper Bail

A long metal bar with two rubber wheels that sit across the reading table. The paper bail holds paper in place during embossing. It is spring loaded for easy removal. Slide the right hand margin to far right, grip paper bail and push to left to remove. It is best removed when a child is writing and reading Braille, giving unimpeded access.

Paper Lever

A yellow plastic lever, located to the right of the embossing area, lifts up and back. When lifted it releases the right margin adjuster and prompts an audio response. When returned to lock position, the audio response prompt says "OK" and the embossing head checks for the size of the page inserted by moving between the margin adjusters.

Ruler

A scale bar at the front of the embossing area to help re-align a Brailled page.

Paper Advance

A yellow plastic grooved roller located to the left of the reading table. Assists with moving a new page into position.

On/Off Switch

Located on the left front is the on/off switch.

Volume Control


At the front of the Mountbatten and to the right of the on/off switch is a thin yellow knob for volume control.

Carry Handle

In the center front is the yellow carry handle.

Pull it out and use it to carry the Mountbatten about.

Location of Parts



Orientation - Ports

Starting from the right-hand side of the Mountbatten:

USE

Volume Control

The USB port on the right side of the case is for making a connection with a PC for Chat, File transfer or embossing.

Carry Handle

Mimic Port

The Mimic port on the right side of the case is for connecting the optional Mimic display.

PC Keyboard Port

On the back right rear of the unit is a round PS2 port for a standard PC keyboard plug.

Parallel-Out Port

Next is a Parallel Out port for sending information to a standard text printer directly from the MB.

Serial Port

The Serial port provides two-way file transfer options to and from a PC or Braille note taker, and MB-Comm software facilitates a simple file transfer program between the Mountbatten and PC.

Speaker

On the left side of the Mountbatten is the speaker for audio output.

Headphone Port

A small socket for a standard headphone jack is located next to the speaker.

AC Adapter Port

An AC Adapter input for mains power supply - this is for the mains charger supplied with every Mountbatten.

DC Adapter Port

A small socket near the corner supports an optional 12 Volt battery charger.



Keyboard Port

Beside the On/Off switch is a port suitable for customised external keyboards to meet special needs.


Location of Ports

Left side profile

Right side profile



Features for loading paper



Embossing head

moves along behind the ruler and embosses Braille dots. The embossing head leads your finger to the ruler.

Paper lever

is on the right - raise it to insert or remove paper.

Reading table

is the flat strip of metal which supports paper.

Metal groove

in the reading table assists to align the top of the page to make sure it is straight

Paper bail

is the round bar that goes from left to the right above the reading table. When Brailleing whole pages the rollers assist with noise reduction. It is spring loaded for easy removal.

Margin adjusters

the left margin is fixed, the right margin adjuster can be adjusted to paper size. Rollers on the margin adjusters grip paper while brailleing.

Manual paper advance

is a knob set in a recess to the left of the paper control area. Use it to adjust the top edge of paper when you are inserting a new page.

Attention!!

Do not use manual paper advance at other times - you will lose the ability to correct or erase.

Loading Paper

- Step One** Lift the yellow paper lever.
- Step Two** Slide the right margin left or right to make room if necessary.
- Step Three** Lay the page flat on the reading table.
- Step Four** Guide the paper under the embossing head.
- Step Five** Push the paper carefully into the Mountbatten.
- Step Six** Check that the paper is placed straight against the left hand margin and under the left grip wheel.
- Step Seven** Align the top of the sheet of paper with the groove at the front edge of the reading table.
- Step Eight** Hold the right side of the paper against the paper roller with your left thumb, and hold the edge up with your left index finger.
- Step Nine** With your right hand, slide the right margin to the left, until it rests against the right edge of the paper. The paper should slide easily under the right grip wheel and over the base of the margin adjuster.
- Step Ten** Flip the yellow paper lever back down. The embossing head will move to check the paper size.

Start-Up

Step One:

Use the information on Ports on Page 3.

Plug in

Plug in the AC power adapter and switch on the power.

If it is not possible to have a power connection during usage, make sure the battery is re-charged after 8 hours of use.

Step Two:

The On/Off switch is on the left front of the Braille.

Switch on

Each time you turn the machine on the embossing head moves to check the paper size

Step Three:

The audio response will greet you and tell you if:

Listen

- the battery is charging
- you are using Learn or Advanced Mode
- if the Mimic is not connected

Note:

If the Mimic is connected audio response will not tell you

Battery Care

The Mountbatten needs power! It uses a lot of power when embossing and has a specially designed battery.


The battery will operate the Mountbatten for up to 8 hours or a day in class. After this ***it must be recharged!***

If possible connect the power adapter provided, at all times. Otherwise, use the power adapter overnight to recharge the batteries.


It is important to remember the following battery care tips:

- ☑ Each time you switch the Mountbatten on, it checks the battery status. It will tell you to charge the battery if it is low.
- ☑ It is a good idea to charge the battery overnight, every night.
- ☑ If the battery is not re-charged when completely run down, permanent damage may result. Replacing a battery is expensive and often inconvenient.
- ☑ It is a good idea to leave the Mountbatten plugged in wherever possible to ensure the battery is kept charged.
- ☑ You need to charge your battery for at least 24 hours prior to storing it over period of up to 6 weeks.
- ☑ To maintain the battery when Mountbatten is out of service for longer periods, call your supplier and arrange to disconnect the battery.
- ☑ You can check battery status by pressing **newline + space + b** - an audio response will tell you the status.

The Mountbatten Keyboard




Braille Keys



- | | |
|---|---|
| ① | ④ |
| ② | ⑤ |
| ③ | ⑥ |

The Swing Cell



The Braille Cell

Using the Mountbatten Keyboard

Moving around the page

tab	Moves the head one tab stop to the right
backtab	Moves the head one tab stop to the left
backspace	Moves the head one space to the left
newline	Moves the head to the left and advances the paper one line
space	Moves the head one space to the right

Erase and Correct

space + backspace	Erase - replaces a character with a space
backspace + the correct dots	Correct - replaces a character with desired character

More moving around the page

newline + backspace	Move one line up the page in the same column
newline + space	Move one line down the page in the same column
newline + enter	Eject the paper

Brailleing on small sized paper


newline + space	Allows Braille writing without paper. The two keys must be pressed while the embossing head is moving.
------------------------	---

Turning Commands On and Off

Learn Mode

In Learn Mode, the command key is disabled in the default setting.

To activate commands




Press the following 4 keys together:

space + newline + tab + backtab

To deactivate press the same 4 keys together. This is called an on/off toggle.

Entering commands

1. Press the **command** key
2. Type the command using Braille keys
3. Press the **enter** key



Writing commands

Commands combine function keys and Braille letters.

For example: entering a command on the Mountbatten, use the Command Key, Braille code, Enter Key; on the PC use the Escape, text and the Enter key.

Why do we use commands?

In Learn Mode

Commands are used in Learn Mode to tell the Mountbatten to:

- set margins, tabs, line spacing and centring
- adjust dot height and paper settings
- turn embossing on and off
- format the Braille page eg page numbers, word wrap
- save files to memory and emboss and delete files
- use graphics mode

In Advanced Mode

Commands are used in Advanced Mode to tell the Mountbatten to:

- continue to perform all Learn Mode functions
- do forward and back translation
- provide structured introduction of Grade 2 Braille contractions [using APH patterns]
- do advanced editing of files in memory
- control synthetic speech parameters

Toggle keys

Command keys are usually 'toggle' keys. This means you repeat the command to undo it.

HINTS

- Never use a space at the beginning or end of a command
- If you make a mistake while entering a command, press the 'command' key to start again.

Speech Commands

For use with recorded and synthetic speech.

spk r	Speaks recorded speech (default) In MB Writer this command indicates the speaking mode is selected
spk a	Speaks all keyboard keys, functions and commands
spk d	Speaks dot keys only – no function keys
spk c	Speaks commands, errors and prompts only – no dot keys
spk n	No speech – tones only
spk s	Speaks synthetic speech (only available in MB Pro)

Formatting Commands

ww	Word wrap
ce	Center text (default off - toggle)
ra	Right adjust (default off - toggle)

Margins

lm	Set left hand margin to 0
lmh	Set left hand margin to embossing head position
lm[n]	Set left hand margin [n = number of columns]
rm	Set right hand margin to 0
rmh	Set right hand margin to embossing head position
rm[n]	Set right hand margin [n = number of columns]
tm[n]	Set Braille top margin [n = number of lines]
bm[n]	Set Braille bottom margin [n = number of lines]

Tab Settings

tsh	Set tab at position of embossing head
tc	Clear all tab settings

Line Spacing

ls 1	Line spacing 1
ls 1.5	Line spacing 1.5
ls 2	Line spacing 2

File Commands

Start a File

`begin [filename]` Opens a new file – filenames can have up to 16 characters

Close a File

`end` Closes an open file

Finding out what's in the memory

`dir` Embosses all file names and their sizes

`free` Embosses all remaining free memory space

Change a File

`Append [filename]` Attaches more information to a file

`rename filename1,
filename2` Renames 'filename1' and calls it 'filename2'

Delete Files

`del [filename]` Deletes a file

`del all` Deletes all files

`undel` Retrieves the last deleted file. This command only works if used before turning the unit off.

Emboss Files from the Memory

`pr [filename]` Prints the file

`mpr n [filename]` Prints 'n' copies of file

Ways to adjust noise level

Embossing sound

The sound of the Mountbatten writing Braille can be minimised in a number of ways:

- place the Mountbatten on a soft surface like 2 computer mouse pads
- use regular photocopy paper and turn impact force or Braille pressure to a minimum
- place the Mountbatten in front of a solid object such as a book or facing wall

Speech sound


At times where a quieter class environment is needed:

- adjust the volume
- use headphones

Adjusting the Braille Impact


For thicker paper and labels

To **increase** the impact, press the Command Key and the Right Function Key together, then release them; one press for every step higher.



For lighter photocopy paper

To **decrease** the impact, press the Command Key and the Left Function Key together, then release them; one press for every step lower.



Each impact level has a number from 0 (no dots) to 9 (highest level) - later Mountbattens will speak this level at each press.

Connecting Mountbatten to a PC with MB-Comm

Switch Mountbatten OFF


Always switch off before making a connection to the MB.

Step by step:

- Load the Mountbatten with paper
- Identify the USB port on the MB and the PC
- Connect the MB and the PC using the cable.
- Find the MB-Comm Icon and run MB-Comm.
- Run MB-Comm and make the connection.


You can make the connection in **one** of two ways:

Click the right-hand red button




OR

On the Device Menu, click Connect Mountbatten



Choosing the Port

A dialog box will appear on the screen:



Important!

If you are using USB


Make sure you choose the correct port before you continue.

The port you choose MUST say "Mountbatten USB".

Click OK when you have chosen the port.

Making sure everything is ready...

Next you will see a reminder:



...and then completing the connection.

If you are happy that everything is ready, click Continue.

What will happen?

The Mountbatten will make a series of short, high-pitched beeps.

The red buttons will turn green one at a time. When the connection is complete, both buttons will be green.

You are now ready to use MB-Comm.

Chat Mode using the PC with MB-Comm

Producing translated Braille and Print simultaneously in real time

Important!

- Mountbatten must be on and loaded with paper
- The PC and Mountbatten must be connected
- MB-Comm must be running and connected
- The two buttons on the MB-Comm screen must be green

Enter Chat Mode

Go to the File Menu in MB-Comm and click "Chat Mode"

- Mountbatten will say: "Back"
This means that the Braille keyboard will send translated text to the MB screen.

Change between users

To change to the PC keyboard, the Braille user must press the **tab** key on the Mountbatten.

- Mountbatten says "Forward"
To change to the Braille keyboard, the Braille user must press the **backtab** key on the Mountbatten.
- Mountbatten says "Back".

Important!

The Braille user on the Mountbatten keyboard

- types first - the Braille keyboard Back Translation is ready first
- changes between Forward and Back Translation using **tab** and **backtab**

The PC keyboard user

- Types text
- produces Grade 2 Braille on the Mountbatten
- produces text on the MB-Comm screen

The Mountbatten keyboard user

- Types Braille
- produces Braille on the Mountbatten
- produces text on the MB-Comm screen

Embossing from the PC with MB-Comm

Producing a translated Braille copy from a text file on the PC

Preparing the file

1. Save the file as Plain Text with a .txt extension
2. Open the file in Notepad
3. Insert any commands you need at the beginning of the file
 - {adv} if the Mountbatten is in Learn Mode
 - {fe} to have Mountbatten translate the file
 - {flit} to preserve print lines in Braille
 - {fgr 1} to have Uncontracted (Grade 1) Braille
 - {m} to emboss at double speed
4. Insert any commands you need at the end of the file
 - {adv} to return to Learn Mode
 - {m} to return to normal speed
5. Save the file again

Important!

Tip:

Insert your commands on a new line by themselves!

Have a folder on the PC to store ready for your text files to be embossed.

Always save to this folder so you can find them again.

See the instructions in previous sections

Make sure you prepare the Mountbatten before you begin!

Run MB-Comm and connect to the PC

From the File Menu select

Emboss file to Mountbatten

Tip

MB-Comm will ask you to choose the file.

Navigate to the file and select it.

Click OK.

Mountbatten will begin embossing the translated file.

Save a Notepad file with only the commands in it to use as a template.

Copy and paste directly from a Word or OpenOffice document into your template.

Embossing Commands

Using Continuous Paper

- `cp` Continuous paper feed (default off - toggle).
- `cpa n` Adds or subtracts 'n' to form length of continuous paper.

Paper Size

- `fl n [.m]` Sets physical length of page in inches or centimetres.

Page Numbering

- `pn [n]` Begins page numbering, starting at 'n'.

Controlling the Braille Output

- `e1 [n]` Allows only 'n' empty lines while embossing to avoid large blank spaces when embossing.
- `cap` Caps lock (default on - toggle).

Using the Mimic

Text display on LCD screen	Print words appear accurately, irrespective of student brailing contracted or uncontracted Braille Text appears after the space bar is pressed
Learn Mode	The Mimic works automatically in Learn Mode. If unplugged or switched off the MB Pro will say "Mimic off". It does not however, say "Mimic on".
Advanced Mode	The Mimic needs to be activated in Advanced Mode.
On/Off switch	The switch is located at the rear left of the unit and needs to be switched on after being connected to the Mountbatten. The Mimic runs on 4 AA batteries or a 6 Volt adaptor.
Automatic Shut down	If not used for 15 minutes the Mimic will shut down. To re-activate, touch any key.
Mimic Keypad	Press keys firmly for more than 1 second to activate. Under the display screen there are three buttons.
Scroll Up	Press left button – one short beep. Scrolls up 20 characters of text.
Scroll Down	Press right button – one short beep. Scrolls down 20 characters of text.
Set Bookmark	Press left and center button – two long beeps.
Clear Memory	Press left and right buttons together for 2 seconds. Will erase the text on the screen.

Forward Translation Commands

Switch to Advanced Mode (command) `adv` (enter)

Turn on Forward Translation

`fe` Forward translate and emboss
`fe [filename]` Forward translate and emboss file
`fm` Forward translate and send to memory
`fm end` End file and turn embossing back on

Change Braille Grade or Code

`fgr [grade]` Change grade [fgr 1 for Grade 1; default Grade 2]
`ffc [language]` Foreign character selection [Refer table next page]

Formatting Your Braille

`flit` Literal mode – a new print line begins a new Braille line
`flit [n]` Creates a hanging indent (n=0 - 5)
`flit s` Suppresses blank lines and adds a 2 cell indent to new paragraphs
`fit` Italics
`fin` Stores current Braille column for indent
`fsh` Multiple space removal

Lists/ Tables

`flist` Creates 2 column table
`flist n1n2` Creates 3 column table

Saving Preferred Format Settings

`save ft` Saves Forward Translation settings
`std ft` Restores all default Forward Translation settings

Changing the Braille Code

ffc Forward translation followed by country code
[language]

bfc Back translation followed by country code
[language]

Country Codes

OZ = Australian
ITL = Italian
DAN = Danish
NOR = Norwegian
UK = English
SP = Spanish
FR = French
SWE = Swedish
GER = German
USA = United States
GRK = Greek
DUT = Dutch
FIN = Finnish
POR = Portuguese

The Options Menu

The Options Menu gives you access to a group of special functions of the Mountbatten Learning System (LS).

To use the Options Menu you must be in Advanced Mode.

Refer to your User Guide for detailed information about the use of each Option.

o-chord

Press the letter o (**dots 1,3,5**) and **space** together.

o-chord is then followed by the appropriate letter to open each of the functions.

z-chord

To exit any of the Options use **z-chord** (**z** plus **space**).

Music

o-chord u

Music Fun mode. Intended for fun and pre-braille.

o-chord u-chord

Music Play Mode. Enter a Braille music character and hear it played back

Stopwatch

o-chord w

Turns on a stopwatch. Test skills, time activities. Have fun.

Time and Date

o-chord d

Have Mountbatten say the current date.

o-chord t

Have Mountbatten say the current time.

o-chord s

Set the current time and date.

Mathematics

o-chord m

Open Simple Math Mode to use the simple calculator.

o-chord n

Change to Nemeth mode. In literary Braille, Nemeth characters will be spoken as numbers.

o-chord c

Turn on a full scientific calculator

The Options Menu - Calculator

Nemeth Mode

o-chord n

Enter Nemeth Mode for practice.

This is NOT a calculator - it changes the speech output so that Nemeth characters are spoken as Nemeth instead of literary braille.

Math Mode or the Simple Calculator

o-chord m

Simple Calculator

- standard calculations including percentage and square root
- interacts with the Editor
- Store and recall and results using memory and variables

Scientific Calculator

o-chord c

Scientific Calculator

- all Simple Calculator features plus complex Scientific and Statistical functions

When you type a calculation you need to use **Nemeth Braille**.

Typing Math into the Calculator - Nemeth Code

Addition	+	⠠	dots 3,4,6
Subtraction	-	⠠	dots 3,6
Multiplication	*	⠠	dots 1,6
Division	/	⠠⠠	dots 4,6 dots 3,4
Percent	%	⠠⠠	dot 4 dots 3,5,6
Square root	√	⠠⠠⠠	dots 3,4,5 number dots 1,2,4,5,6
Decimal point	.	⠠	dots 46
Left parenthesis	(⠠	dots 1,2,3,5,6
Right parenthesis)	⠠	dots 2,3,4,5,6
Comma	,	⠠	Dot 6

Simple fraction

dots 1,4,5,6 number dots 3,4 number dots 3,4,5,6

:::number ::: number ::: – with no spaces between elements.

Mixed fraction

dots 3,4,5,6 number dots 4,5,6 dots 1,4,5,6 number dots 3,4 number dots 4,5,6 dots 3,4,5,6

:::number ::: :::number ::: number ::: ::: – no space between elements.

Numbers in Nemeth Mode

1	.
2	::
3	::.
4	:::
5	::.
6	:::
7	:::
8	::.
9	::.
0	:::

Calculator Commands

e-chord

Execute the calculation or speak the current result.

c-chord

Speak the current line.

t-chord

Emboss the current equation.

dots 3,5,6-chord

Clear the Calculator.

or

zero-chord

f-chord

Speak numbers as digits or words. - this is a toggle.

p-chord followed by a digit from 0 - 6, then	Set the precision level - the number of significant digits past the decimal point used in a calculation.
e-chord	
s-chord followed by any letter from a to z except r	Store the current result in memory at a location with the name of the letter chosen.
dots 3,4,6-chord also known as plus-chord	Change the sign of the result to positive.
dots 3,6-chord also known as minus-chord	Change the sign of the result to negative.
o-chord k	Calculates the expression on the current line of the currently open file in the Editor.
o-chord i	Insert the result of the current calculation into the Editor.
d-chord	Set the argument of commands as degrees
r-chord	Set the argument of commands as radians

Calculations in the Editor

In the Editor

Go to a new, empty line in your document

Write your calculation all on the one line, without pressing newline

- Make sure you use all the correct Nemeth operators
- You can use spaces if you want

Press **o-chord k** to invoke the Calculator and hear the result

Press **o-chord i** to insert the result into your document

The Options Menu - Music


Music – Simple Mode

o-chord u Simple Music mode. Intended for fun and pre-braille.

The six dot keys and the command key become a scale. Press each key sequentially:

3 2 1 command 4 5 6

Press the same keys again to hear the scale in a higher octave.



- **Change Pitch**
- **Sharpen or Flatten**
- **Notes of the Scale**
- **Change Duration**

Other keys have the following functions:

backspace	Move scale down an octave
enter	Move scale up an octave
space	Extend the note by a quarter bar/crotchet
newline	Extend the note by a half bar/minim
newline + space	Extend the note by a full bar/semibreve
tab	Flat
backtab	Sharp

Press multiple note keys together to hear a chord, maximum three notes.

z-chord Exit Music Fun Mode.

o-chord u-chord Music Play Mode.

Enter a Braille music character and hear it played back

In Music Play mode, students can type Braille music characters and hear them being played at the same time. This mode is intended for students who are just starting out learning Braille and becoming familiar with the Braille music code.

Braille Music Notes

C	D	E	F	G	A	B	Rests	Type
⠠	⠡	⠢	⠣	⠤	⠥	⠦	⠇	Wholes or 16ths
⠧	⠨	⠩	⠪	⠫	⠬	⠭	⠮	Halves or 32nds
⠏	⠐	⠑	⠒	⠓	⠔	⠕	⠖	Quarters or 64ths
⠗	⠘	⠙	⠚	⠛	⠜	⠝	⠞	8ths or 128ths


Braille Graphics

Making Braille Graphics

g£


Reduces the space between dots so that lines can be drawn.

Type these simple Braille graphics from Sheila Amato on the Mountbatten keyboard.




Heart

1. s, e, i, wh
2. gh, space twice, ar
3. space once, e, i



Teddy Bear

1. space once, ow, o, two colons, ow, o
2. space once, s, ow, o, ow, o, wh,
3. space once, gh, space once, e, i, space once, ar
4. ow, s, space 4 times, wh, o
5. in, gh, space 4 times, ar, en
6. e, i, e, 2 colons, i, e, i



Sailboat

1. space once, 4-5-6, ch
2. space once, 4-5-6, space once, ch
3. space once, 4-5-6, space twice, ch
4. er, full cell, full cell, full cell, full cell, n

Writing Numbers in a Command

Option 1 Use the number sign, then the number

Option 2 Use the dropped Nemeth number

Braille Writing Commands

The following commands apply in both Learn Mode and Advanced Mode

cor	Correction Mode (When off, correcting and erasing are disabled)
em	Turn embossing off (default on - toggle)
man	Choose between manual and automatic new line (default auto - toggle)
ms [n]	Multi strike feature [n=1-4] can be used to Braille heavy materials
np	New Page or page eject
pn [n]	Begin page numbering
pbell [n]	Set page bottom warning bell to ring [n] lines before bottom of page
rmbell	Right Margin Bell (rings 5 positions before end of line) (default on – toggle)
show	Emboss commands as they are entered (default off – toggle)

Miscellaneous Commands

newline + space + b	Checks the battery status
ku	One Handed Mode. The default status is off; ku is a toggle command.
ver	Embosses or speaks the Mountbatten software version.
hp	Help. Embosses a command summary.
save	Saves all current settings
std	Restores the Mountbatten factory settings
alt n	Enters non-standard ASCII characters

Making Braille Graphics

gf	Reduces space between dots so that lines can be drawn. The default status of gf is off; this is a toggle command.
-----------	--

Adjusting Dot Height

backtab + command	Increases the embossing impact
tab + command	Decreases the embossing impact
ms [n]	Multi-strike. (n=1 to 4)

Opening the Serial Port

s1	Opens the Serial Port so that text can be read on the Mimic while a PC keyboard is in use [forward translation]
-----------	---

Braille characters and contractions

This table shows the BANA Braille characters and contractions used in North America.

a	b	c	d	e	f	g	h	i	j
k	l	m	n	o	p	q	r	s	t
u	v	x	y	z	and	for	of	the	with
ch	gh	sh	th	wh	ed	er	ou	ow	w
/ea/, '	be/ bb/ ;	con/ cc / :	dis/ dd / .	en	to / ff / !	were (/gg/)	his "/>	in	was by/ /
st	ing	num/ ble/ ble	ar	'	com -				
dot 4	dot 4,5	dot 4,5,6	dot 5	dot 4,6	letter sign	cap			

Student Evaluation Tool

For Teachers and Para-professionals

This evaluation tool gives educators a basic skills checklist that can help in evaluating skills developed as the student learns to use the Mountbatten for every-day activities.

The list can be used to set achievable objectives when drawing up instruction plans. Dating when skills were acquired will assist in measuring growth and provide a basis for review.

Space has been left at the bottom of each section for comments to assist with monitoring.

Student Name:	Class:
School:	
Assessing Teacher:	Date:

1. Mountbatten Orientation – Parts:

OBJECTIVE	YES	DATE
Student can turn Mountbatten on and off		
Student knows location of all 6 Braille keys		
Student knows the dot number for each key		
Student can operate volume control		
Student can locate the following parts of the Mountbatten:		
Paper lever		
Margin adjusters [left and right]		
Embossing head		
Paper advance knob		
Student knows the location of the blue keys:		
Space bar		
Back space		
Newline		
Enter key		
Command		
Tab		
Backtab		
Student knows how to care for equipment:		
No food or drink near equipment		
Clean hands before use		
Turn off after use		
Use of power supply		
Care of battery		
Keep away from desk edge		
Location of carry handle		
Portability care		

Comments:

2. Mountbatten Orientation – Ports

OBJECTIVE	YES	DATE
Student can locate the following ports and identify the appropriate plug: Headphone port and jack		
Power adapter port and plug		
PC keyboard port and plug		
Serial and Parallel printer ports		
Serial cable and plug		
Parallel cable plug		
Mimic port		
Mimic cable and plugs		
Student will understand the concepts:		
Left front		
Right front		
Left back		
Right back		
Left side		
Right side		

Comments:

3. Set-up:

Note: Being able to successfully plug in all connecting devices is a fine motor skill which requires development through assistance and instruction. This set of Objectives along with Objective 8, 9 and 10 may best be assessed when the skills required are developed rather than when they are first introduced.

OBJECTIVE	YES	DATE
Student can connect the power adaptor to the Mountbatten and		

switch it on		
Student can switch Mountbatten on and off		
Student is aware of content of recorded speech output on Mountbatten at start-up (with and without Mimic connected)		
Student is aware of recorded speech response when Brailleing		
Student can connect the headphones and use the volume control		

Comments:

4. Loading Paper

OBJECTIVE	YES	DATE
Student can remove and replace paper bail		
Student is aware of audio response when paper lever is lifted		
Student can load photocopy paper from the top		
Student can load Braille paper from the top		
Student can align paper with groove on reading table and position it correctly under wheel of margin adjusters		
Student can decrease Braille pressure		
Student can increase Braille pressure		
Student can remove paper manually		
Student can use automatic page eject		
Student can correctly load a pre-brailled worksheet		

Comments:

5. Brailleing in Grade 1

OBJECTIVE	YES	DATE
Student can place fingers on Braille keys correctly		

Student can use space key and new line key		
Student can type Grade 1 Braille		
Student can use audio response to support action		
Student can identify and understand when to use Command key		
Student can use reading table to track progress		

Comments:

6. Moving around the page:

OBJECTIVE	YES	DATE
Student can erase		
Student can erase and correct		
Student can move down page without moving embossing head		
Student can move up the page without moving embossing head		
Student can move across the page and back to left margin		
Student understands the use of and can switch on Word Wrap		
Student can centre a word		
Student can set double line spacing		
Student can widen left margin		
Student can insert page numbers		
Student can execute a command sequence accurately		

Comments:

7. Creating and Saving Files in Memory

OBJECTIVE	YES	DATE
Student can start and name a file		
Student can end a file		

Student can emboss a directory		
Student can select a file name		
Student can emboss the contents of a stored file		

Comments:

8. Connecting Devices – Mimic

Note: *Being able to successfully plug in all connecting devices is a fine motor skill which requires development through assistance and instruction. This set of Objectives may best be assessed when the skills required are developed rather than when they are first introduced.*

OBJECTIVE	YES	DATE
Student can identify Mimic port and connect the Mimic		
Student can switch on Mimic		
Student is aware of changes to audio output on start-up		
Student knows to press space to send information to Mimic		
Student knows that in Learn Mode text is displayed automatically		
Student is aware of the automatic shut-off feature		
Student is aware of and can activate the wake-up feature		
Student is familiar with the functionality of the Mimic		

Comments:

9. Connecting Devices – PC keyboard

Note: *(1) Being able to successfully plug in all connecting devices is a fine motor skill which requires development through assistance and instruction. This set of Objectives may best be assessed when the skills required are developed rather than when they are first introduced.*

(2) Learning the layout of the QWERTY keyboard is a separate skill development exercise.

OBJECTIVE	YES	DATE
Student is aware of the correct procedure when connecting a PC keyboard [switch off Mountbatten prior to connecting PC keyboard]		
Student can identify the correct port and plug in the PC keyboard		
Student is aware of the main features of Advanced Mode		
Student can activate Advanced Mode and return to Learn Mode		
Student understands the concept of Forward Translation		
Student knows the command for Forward Translation		
Student can listen effectively to synthetic speech		
Student is aware that Grade 2 Braille is the automatic default in Advanced Mode		
Student is aware of the advantages of the PC keyboard for sighted users		

Comments:

10. Connecting Devices – Printer

Note: Being able to successfully plug in all connecting devices is a fine motor skill which requires development through assistance and instruction. This set of Objectives may best be assessed when the skills required are developed rather than when they are first introduced.

OBJECTIVE	YES	DATE
Student is aware of correct procedure when connecting a printer [Switch off Mountbatten prior to connecting printer]		
Student can identify the correct port and can plug in the printer		
Student can identify plugs and insert them correctly		
Student can check for paper in printer		
Student understands the concept of Back Translation		
Student can turn on Back Translation		
Student can send a file to be embossed & printed simultaneously		

