

Important note: All technical specifications are subject to change without notice.

Unit Dimensions and Weight

Item	Height	Width	Depth	Weight
explorē 8	0.94"	8"	6"	1.30 lbs – 20.8 oz
	24 mm	204 mm	153.2 mm	580 g
explorē 8 folding stand (folded)	1.6"	13.6"	10.6"	2.44 lbs – 39.04 oz
	39 mm	345 mm	270 mm	1105 g
explorē 8 folding stand (opened)	10.3"	13.6"	10.6"	2.44 lbs – 39.04 oz
	261 mm	345 mm	270 mm	1105 g

Packaging Dimensions and Weight

Item	Height	Width	Depth	Weight
explorē 8 (device)	3.15"	9.76"	7.09"	2.70 lbs – 43.2 oz
	80mm	248mm	180mm	907 g
2.7explorē 8 Kit (device, stand, bag)	8.5"	16"	13.5"	8.5 lbs – 136 oz
	210 mm	410mm	343mm	3856 g
explorē 8 Stand (stand, bag)	7"	16"	13.5"	5.5 lbs – 88 oz
	178 mm	406mm	343 mm	2495 g

explorē 8 Specifications

Specifications	explorē 8
Screen	8-inch IPS LCD screen
Battery	Non-removable Li-Polymer 6700 mAh battery
Camera	Two HD 21M pixels
Power Input	C-Type USB connector, 5V@3A and Quick charge 12V@1.5A
Charge Time	2h
HDMI Compatibility	Supported HDMI formats: 1280x720p 60Hz, 50Hz 1920x1080p 60Hz, 50Hz, 24Hz 1920x1080i 60Hz et 50Hz No support of DVI emulation.

explorē 8 Features

Features	explorē 8
Zoom	2x – 30x
Focus	Close-up: Fix auto-focus Far: auto-focus or locked auto-focus Stand ONLY: Temporary focus lock
Languages available	English US (is default), English UK, French, Spanish, German, Dutch, Danish Polish, Portuguese (Portugal), Portuguese (Brazil), Norwegian(Bokmål), Italian, Swedish, Finish, Russian, Czech, Japanese, Korean, Icelandic, Estonian, Turkish, Lithuanian, Hungarian, Latvian, Arabic
Touch screen gesture configuration	Live zoom gestures, tap & hold to save
Preferred zoom	Save a preferred zoom level
Line and blinds	Enable/disable line or blinds
Menu text size	Small, Medium and large
Battery Life	4-hours
Power saving	Auto shutoff after 1min, 2mins, 5mins or 10mins, never.
Color schemes	Full color, White on black, Black on white, Yellow on black, Black on yellow, Orange on black, Black on orange, Black on cyan, Cyan on black, Black on purple, Purple on black, Black on green, Green on black, Blue on yellow, Yellow on blue, White on blue, Blue on white, Red on white, White on red, Red on black, Black on red
Display frequency	50Hz-60Hz
Close-up viewing distance	8.5 cm – 15 cm / 3.35” – 5.90”
Far viewing distance	15cm – infinity / 5.90” - infinity

explorē 8 Operating Conditions

Specifications	explorē 8
Temperature	10 °C to 40 °C / 50 °F to 104 °F
Humidity	20% to 80%

explorē 8 Box Contents

explorē 8 Device Contents
explorē 8 electronic magnifier
USB Charger
USB Cable
HDMI cable
User manual
Quick Start guide
Warranty Card
Product Registration Card
Protective case
LCD Protector
explorē 8 Kit Adds
Stand
Carry Bag

explorē 8 Stand Accessory Contents
Stand
Carry Bag
Quick Start Guide
Stand Notice