



Supported Optical Character Recognition and Text-to-Speech

| OCR 52 Languages | TTS 16 Languages | TTS Details 77 Voices; 16 Dialects |
|----------------------------|----------------------------|---|
| Afrikaans | | |
| Albanian | | |
| Basque | | |
| Belarusian-Cyrillic | | |
| Breton | | |
| Bulgarian-Cyrillic | | |
| Corsican | | |
| Croatian | | |
| Czech | ✓ | 1 |
| Danish | ✓ | 3 |
| Dutch | ✓ | Belgium-4, Netherland-4 |
| English | ✓ | US-10, Australia-2, Canada-1, India-2, Seychelles-1, UK-5 |
| Estonian | | |
| Faroese | ✓ | 2 |
| Finnish | ✓ | 1 |
| French | ✓ | Canada-3, France-8 |
| Frisian | | |
| Gaelic-Irish | | |
| German | ✓ | 5 |
| Greek | ✓ | 1 |
| Hungarian | | |
| Icelandic | | |
| Italian | ✓ | 3 |
| Kurdish | | |
| Latvian | | |
| Lithuanian | | |
| Luxembourgian | | |
| Maltese | | |
| Moldavian-Cyrillic | | |

| | | |
|--------------------|---|----------------------|
| Norwegian | ✓ | 3 |
| Polish | ✓ | 2 |
| Portuguese | ✓ | Brazil-2, Portugal-1 |
| Provençal | | |
| Romanian | | |
| Sami | | |
| Sami-Lule | | |
| Sami-Northern | | |
| Sami-Southern | | |
| Serbian-Cyrillic | | |
| Servian-Latin | | |
| Slovak | | |
| Slovenian | | |
| Sotho | | |
| Spanish | ✓ | Spain-4, US-2 |
| Swahili | | |
| Swedish | ✓ | Sweden-5, Suomi-1 |
| Tahitian | | |
| Turkish | ✓ | 1 |
| Ukrainian-Cyrillic | | |
| Welsh | | |
| Xhosa | | |
| Zulu | | |
| NOTE: | <p>Prodigi supports OCR in many languages, but not every OCR language has a matching Acapela text-to-speech (TTS) voice. When no Acapela voice is available, you may install a 3rd party SAPI-5 TTS voice for that language. These SAPI-5 voices can read document files only.</p> <p>Open a file → Settings → Reading Voice to select an installed SAPI-5 TTS voice.</p> <p>System menus, messages, and the user interface will always be read using Acapela voices.</p> | |