



Prodigi software

powered by Humanware

Prodigi for Windows

***Technical Requirements and Specifications**

Software Requirements: The following are suggested *minimum requirements* for using Prodigi software with user-supplied hardware.

Computing Device (tablet, laptop, or PC)

Description	Specification
Operating System (64-bit)	Windows 11, Windows 10
Processor*	Intel Core i3 (4 CPUs) or better CPU
Memory (RAM)	At least 8 GB
Required hard disk space	8 GB free disk space
Video	DirectX 11 or newer
Sound	Windows-compatible sound card (for speech).

*ARM64 may not support the optional Kodak distance camera.

Document Camera

Description	Specification
Resolution	13MP
Color	24-bit color
FPS	60fps@10ADP
Operating System	Windows 10+
Connectivity	USB connectivity

Internet: Required for installation, activation, some features (example – Books), and ongoing updates.

URL Whitelist: The following URLs must be whitelisted for proper functionality in restricted networks:

- <https://hgphcardsad01.blob.core.windows.net/>
- <https://id.humanware.com>
- <https://humanware.ca/>



Prodigi software

powered by HumanWare

****Hardware Specifications:** The following are the basic specifications for the hardware sold with the Prodigi Starter Kit, Complete Kit, and All-in-One configurations.

Hovercam Solo 8+

Description	Specification
Resolution	4224 x 3156; 13MP
Color	24-bit color
FPS	4K @30fps, 1080p @60fps, 720p @120fps
Operating System	Windows 10+
Connectivity	USB 3.0
Dimensions	1.5 x 3.5 x 11.4" 3.8 cm x 9 cm x 29 cm Extended: 19.56" /49.75 cm
Weight	2.2 lbs /1.0 kg
Service & Support	Hovercam 2-year parts and labor warranty

[Click here for complete Solo 8+ specifications.](#)

Microsoft Surface Pro 9

Description	Specification
Operating System	Windows 11 Pro
CPU	Intel Core i5 (12th Gen) 1245U / 1.6 GHz
Memory (RAM)	8GB; LPDDR5 SDRAM
Storage	128GB SSD
Display	13" touchscreen (10-point multi-touch) 2880 x 1920
Audio & Video	Intel Iris Xe Graphics Cameras: 4K (rear); 1080p (front) Resolution: 10 Megapixel (rear) Stereo speakers, dual far-field studio microphones Dolby Atmos compliant
Wireless	Bluetooth 5.1, 802.11a/b/g/n/ac/ax (Wi-Fi 6E)
Battery	Up to 15.5 hours
Connections	2 x USB4 (supports Thunderbolt 4) (Type C)
Dimensions	WxHxD: 11.3 x 8.2 x 0.4" WxHxD: 28.7 x 20.8 x 1.0 cm
Weight	31oz. / 878.8g
Service & Support	Microsoft Limited warranty - 1 year

[Click here for full Microsoft Surface Pro specifications.](#)



Prodigi software

powered by HumanWare

AOC 24" All-in-One PC (available in North America only)

Description	Specification
Operating System	Windows Pro 11
Processor	Intel® Core™ i3-12100
Memory	DDR 16GB
Storage	512GB SSD
LED Monitor	23.8", IPS Full HD, 1920x1080, matte
Network	Realtek 8111EL&F 100/1000M Ethernet, RJ45 port
Audio	Realtek ALC662, Azalia Code support 5.1 channel output, 2 Jacks, MIC in/Line out
Wi-Fi	802.11b/g/n/ac 2.4GHz & 5GHz Dual Band WIFI 5 and Bluetooth 4.0
Speakers	Two 3W Height adjustable stand with quick release function
Webcam	3MP
Keyboard and Mouse	AOC USB
Power Adapter	13.5/19V
UPS Battery	60 minutes
Input/Output	1 line out and 1 Mic in, COM Port, 2 - USB 2.0 Ports, 4 USB 3.0 Ports, LAN RJ45, HDMI, DC-In port
Dimensions	Monitor: 23.9" diagonal and (HxW) 13.5 x 21.25" Monitor: 60.7cm diagonal and (HxW) 34.3 x 54.0cm Assembled PC (DxHxW): 8 x 14.5-19.5 x 21.25" Assembled PC (DxHxW): 20.3 x 36.8 x 54cm
Weight	Assembled PC: ~11.9lbs / 5.4kg Assembled PC, keyboard, and mouse: ~18.0lbs / ~8.2kg

**Specifications and requirements are subject to change at any time.*

***HumanWare reserves the right to substitute hardware with equal or greater specifications at any time.*