

*Important note: All technical specifications are subject to change without notice.*

### Unit Dimensions and Weight

Item	Height	Width	Depth	Weight
Reveal16 Open	14.57"	14.96"	18.3"	13 lbs
	370 mm	380 mm	464 mm	5900 g
Reveal16 Folded	5.9"	14.96"	18.6"	13 lbs
	150 mm	380 mm	473 mm	5900 g
Reveal 16 with X-Y Folded	7.1"	14.96"	18.6"	Adds 5.74 lbs
	180 mm	380 mm	464 mm	Adds 2604 g
Optional Battery				Adds 0.89 lbs
				Adds 404 g

### Reveal 16 Specifications

Screen Specifications	Reveal 16
Display	15.6" diagonal LCD
Aspect Ratio	16-9
Resolution	1920x1080
Minimum Viewing angle	70
Pixel Arrangement	RGB Vertical Strip
Power Consumption	3.08W
Display frequency	50Hz, 60Hz

Camera Specifications	Reveal 16
Resolution	1080p
Minimum Live Magnification	1x
Maximum Live Magnification	45x
Minimum Optical Magnification	1x
Maximum Optical Edge Magnification	10x
Refresh Rate	60 fps
Working Distance Between Work Surface and Camera Housing	11.8" / 30 cm
Operational distance for distance viewing	100 cm to infinity
Focus	Auto-focus

Battery Specifications	Reveal 16
Battery autonomy	5 to 6 hours with Lighting on.
Battery charging time	6 hours
Maximum Charge/discharge Current	6800mA
Battery Indicator	Yes

### Reveal 16 Features

Features	Reveal 16
Buttons	4: Power, Auto-Focus Lock, Zoom, Contrast Colors
Power save/sleep Mode	Yes
On/Off Switch	Yes
Brightness Control	Yes
Magnification Control	Yes
True Color	Yes
Color Select mode	Yes
Enhanced Contrast	Yes
Auto-Focus Lock	Yes
Manual Auto-Focus	Yes
Document Overview	Yes
Camera and Color Contrast Threshold adjustments	Yes
Adjustable Line and blinds	Yes
Enhanced Contrast Options	Full Color, Programmed Contrast Colors. (Default: white on black, black on white, yellow on black and black on yellow)
Color schemes	Full color, Black on white, White on black, Black on yellow, Yellow on black, Black on orange, Orange on black, Black on cyan, Cyan on black, Black on purple, Purple on black, Black on green, Green on black, Blue on yellow, Yellow on blue, White on blue, Blue on white
Menu Languages	English US (is default), English UK, French, Spanish, German, Dutch, Danish Polish, Portuguese (Portugal), Portuguese (Brazil), Norwegian (Bokmål), Italian, Swedish, Finish, Russian, Czech, Japanese, Korean, Icelandic, Estonian, Turkish, Lithuanian, Hungarian, Latvian, Arabic
HDMI Compatibility	Supported HDMI formats: 1280x720p 60Hz, 50Hz 1920x1080p 60Hz, 50Hz, 24Hz 1920x1080i 60Hz et 50Hz No support of DVI emulation.

### Reveal 16 Box Content

Content
Reveal 16 Device
User Guide
Getting Started Sheet
Power Adapter
North American Power Cable
European Power Cable
UK Power Cable
Australia Power Cable